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Manufacturer: Beijing Safe Heart Technology Ltd.

Address: Room 101, Unit 6, Building No.6, No.88 kechuang 6th Street, Beijing Economic-Technological Development Area, 101111 Beijing, P.R. China

European Representative:
Shanghai International Holding Corp. GmbH(Europe)
Address: Eiffeistrasse 80, 20537 Hamburg, Germany
Email: shholding@hotmail.com

General Description


Oxygen saturation (SpO₂) is a percentage of oxygen contained in the blood. 'S' indicates saturation, 'p' indicates the pulse and O₂ is the symbol for oxygen. Blood cells contain a chemical called haemoglobin (Hb). When haemoglobin is combined with oxygen it forms a second chemical called oxyhaemoglobin (HbO₂). As blood travels around the body the oxyhaemoglobin releases the oxygen where required and then turns back into haemoglobin. Oxygen levels in the blood can fluctuate, especially during activity. Generally normal oxygen levels of a person is 95% or greater. Less than 95% may indicate health issues which would require clinical investigation. Oxygen deficiency in the blood is referred to as hypoxia. There are many respiratory conditions that can result in a patient becoming hypoxic; COPD (Chronic Obstructive Pulmonary Disease), OSA (Obstructive Sleep Apnea) and chronic bronchitis to name a few. Oxygen saturation plays an important part for clinical professionals to be proactive in patient diagnosis.

The Drive DeVilbiss Healthcare Hb0-2000 fingertip pulse oximeter is compact in size, user friendly and features an OLED screen display. SpO₂ diagnosis requires the patient to simply wear the device on their finger. The oxygen saturation and pulse rate values will be displayed on the device screen. The Drive DeVilbiss Healthcare Hb0-2000 has been designed and tested to provide repetitive, accurate screening for SpO₂ patient diagnosis.

Measurement principle

The principle measurement of an oximeter uses a mathematical formula known as the 'Lambert Beer Law' according to spectrum absorption characteristics of reductive hemoglobin (RHb) and Oxyhemoglobin (HbO₂) in glow and near-infrared zones. Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with capacity pulse scanning and recording technology, so that two beams of lights at different wavelengths (660nm glow and 905nm near infrared light) can be applied to the fingertip via a clamp finger-type sensor. A measured signal obtained by a photosensitive element, and processed via the device electronics and microprocessor will be displayed on the oximeter's screen display.

- Diagram of Operation Principle
1. Red and Infrared-ray Emission Tube
 2. Red and Infrared-ray Receipt Tube

Product Operation Scope

The Hb0-2000 fingertip pulse oximeter is a non-invasive device designed for oxygen saturation (SpO₂) and pulse rate detection of adult and pediatric patients in the home and clinical environments. It has not been designed for continuous monitoring.

Precautions for use

1. Carefully read this user manual before using the Hb0-2000 fingertip pulse oximeter
2. Do not use the Hb0-2000 fingertip pulse oximeter in an MRI or CT environment.
3. Do not use the Hb0-2000 fingertip pulse oximeter in situations where alarms are required. The device has no alarms.
4. A pulse oximeter is intended only as a supplementary accessory in patient assessment. It must be used in conjunction with other methods of clinical diagnosis.
5. Frequently check application site to determine correct positioning of the sensor.
6. Check patients circulation before use and refrain from using if the skin becomes sensitive.
7. The pulse oximeter has no alarms; it is not intended for continuous monitoring.
9. Prolonged use of the device or the patient's condition may require changing the site of the sensor periodically.
10. Autoclaving, ethylene oxide sterilizing, or immersing the Hb0-2000

fingertip pulse oximeter in liquids will damage the device

11. The following conditions could cause inaccurate readings:


- High levels of dysfunctional haemoglobins such as carboxyhaemoglobin or methemoglobin
- Intravascular dyes such as indocyanine green or methylene blue.
- High ambient light.
- Excessive patient movement.
- High-frequency electrosurgical interference.
- Venous insufficiency.
- Application of a blood pressure cuff, arterial catheter, or intravascular line.
- Hypotension, severe vasoconstriction, severe anemia, or hypothermia.
- Cardiac arrest or shock.
- Fingernail polish or cosmetic fingernails.

Please comply with local authority regulations when disposing/recycling of the device and batteries.

Product features

1. The Hb0-2000 Fingertip pulse oximeter is light weight, compact and easy to use.
3. Low power consumption only requires 2 X AAA batteries. (20 hours)
4. A low voltage warning will be displayed on the device screen when batteries are low.
5. The product will automatically power off after 8 seconds if no signal is detected.
6. No maintenance required.

Operation Instructions

The Hb0-2000 features six different directional screen settings including a pulse rate bar and plethysmography waveform.

Parameter setting:

Press the power button (> 0.5s), the oximeter will enter into the parameter setting screen. There are two parameter sub menu screens:

Devices with sound/alarm function

When the "*" signal is shown on the "Sounds Setup", press the button (> 0.5s) and enter into the sounds setting menu, pressing the button in turn will select the item to be changed, Press and hold the button to change the parameter. Select "+" or "-" to increase/decrease the number values.

When the "*" signal is shown on the "Alm Setup", press the button (> 0.5s) and enter into the alarm setting menu. Press the button to set the on/off for the alarm and beep.

When the "*" signal is shown on the "Restore", press the button (> 0.5s) to default to factory settings.

Devices without sound/alarm function

SpO₂ and PR values can still be set for units without sound. When the SpO₂ and PR increase/decrease outside of the parameter settings the measured unit will flash to indicate a warning.


Settings	Settings
Alm Setup *	Sounds Setup *
Alm off	SpO ₂ Alm Hi 100
Deep off	SpO ₂ Alm Lo 90
Restore OK	PR Alm Hi 130
Exit	PR Alm Lo 50
	Exit

Battery Installation

1. Open the battery compartment by sliding the battery door
2. Place the two AAA batteries into the battery compartment ensuring correct position.
2. Slide the battery door back in position ensuring battery compartment is fully closed.

Important:

- To prevent damage ensure correct polarity when installing/replacing batteries.
- To prevent damage install negative (-) side of battery first.
- For extended periods of non-use remove the batteries from the oximeter

Lanyard Installation

1. Thread the thin end of the lanyard through the hole at the rear of the case.
2. Loop the wider end of the lanyard through the thin end before pulling it tight.

Operation Instructions

1. Install two AAA batteries (See battery installation)
2. Insert finger into the oximeter as shown in Pic
3. The finger nail must be facing upwards.

4. Switch the device on by pressing the power button.
 5. Avoid excessive body movement during measurement.
 6. Measurement data will be displayed on the device screen.
 7. Remove finger to stop measurement. Device will automatically switch off after 8 seconds.
- Noximeter between patients. The Hb0-2000 oximeter is constructed of latex free materials.

Product Specifications

1. Display Type: OLED (Organic Light Emitting Diode)
2. SpO₂:
Measurement range: 35%-100%
Accuracy: 70% - 99% $\pm 3\%$ $\leq 69\%$, no definition
3. Pulse Rate:
Measurement range: 30-250BPM
Accuracy: 30 - 250BPM $\pm 2\%$
Pulse rate display model: Bargraph Indicator
4. Power Requirements:
Two AAA alkaline batteries
Power consumption: Less than 30mA
Low power indication: The low voltage warning will illuminate when the voltage is less than 2.6V($\pm 0.2V$)
Battery Life: Two AAA 1.5V, 600mAh alkaline batteries 20 hours operation.
5. Dimension:
LxWxH: 57mm x 31.5mm x 30.5mm
6. Environment Requirements:
Operation Temperature: 5 - 40°C
Storage Temperature: 20- 70°C
Ambient Humidity:
work: 15%-90%; Storage: 15%-93%
7. EMC of this product complies with IEC60601-1-2 standards.
8. Measurement performance in low perfusion condition: Required test equipment (FLUKE-INDEX 2XL) the pulse wave is functional without failure when the simulation pulse wave amplitude is set at 0.4%.
9. Interference Resistance Capacity against ambient light: Device works normally. Mixed noise produced by FLUKE-INDEX 2XL

Product Accessories

1. Lanyard
2. 2 X AAA batteries
3. User manual

Maintenance and Storage

- Replace the 2 X AAA batteries when low voltage warning is indicated.
- Clean the internal/external surfaces of the Hb0-2000 fingertip oximeter before and after patient use.
- For extended periods of non-use remove the batteries from the oximeter
- Storage conditions: Ambient temperature 10°C -40°C (14-104°F) and relative humidity 10%-80%
- Store the unit in a damp free, dry environment.
- Please comply with local authority regulations when disposing/recycling of the device and batteries.

Calibrating the pulse oximeter

Factory calibration is conducted using a FLUKE-INDEX 2XL oxygen analogue meter, and BCI oximeter calibration meter.

Note: The test methods used to establish the SpO₂ accuracy requires

clinical testing. The oximeter used to measure the arterial haemoglobin oxygen saturation levels and these levels are to be compared to the levels determined from arterial blood sampling with a CO-oximeter. For more information on device calibration please contact your local Drive DeVilbiss Healthcare representative.

Troubleshooting

Problems	Possible reason	Solution
SpO ₂ or Pulse Rate not being displayed	Finger is not inserted correctly	Re-insert or try a different finger
SpO ₂ or Pulse Rate not being displayed	Excessive body movement	Reduce and stabilise body movement
Reduce and stabilise body movement	Exhausted batteries Batteries incorrectly installed Faulty oximeter	Replace batteries Reinstall the batteries Please contact your local customer service centre
Device display turns off suddenly	The product will automatically power off after 8 seconds if no signal is detected Exhausted batteries	Normal Replace the batteries

Symbol Definitions

	Type B applied part
	Attention, refer to the manual before using
%SpO ₂	Oxygen saturation
PR _{bpm}	Heart rate (BPM)
	Low power indicator
	No SpO ₂ alarm
SN	Serial number

Manufacturer: Beijing Safe Heart Technology Ltd.

Address: Room 101, Unit 6, Building No.6, No.88 Kechuang 6th Street, Beijing Economic-Technological Development Area, 101111 Beijing, P.R. China

European Representative:
Shanghai International Holding Corp. GmbH(Europe)
Address: Eiffelstrasse 80, 20537 Hamburg, Germany
Email: shholding@hotmail.com

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Manufacturer:Beijing Safe Heart Technology Ltd.