Prüfanweisungen Badelifter

BELLAVITA, FELIX UND FELIX T

Zweck

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BELLAVITA, FELIX UND FELIX T

Intention

According to MPBetreibV operators of medical devices are bound to perform safety-related tests, that follow generally accepted engineering rules. To determine failures, which might occur once in a while according to experience, in time, safety-related tests have to be carried out in adequate periods. We recommend to do this every 12 months.

Product description

Product Name: _____

Serial No.: _____

Date of manufacture: _____

Customer: _____

Premises

Safety-related test must be carried out by a trained and experienced person, who can guarantee a correct carrying out of the test, who isn't subject to any directive according to that control task and who has measuring and testing devices at her/his disposal. An official documentary is required for every test. It must be retained until the next safety-related test take place.

General condition inspection

- | Ja | Nein | |
|--------------------------|--------------------------|--|
| <input type="checkbox"/> | <input type="checkbox"/> | Is the product complete with all parts? |
| <input type="checkbox"/> | <input type="checkbox"/> | Has there been any changes at the product? |
| <input type="checkbox"/> | <input type="checkbox"/> | Is the product free of damages? |
| <input type="checkbox"/> | <input type="checkbox"/> | Is the intended function given? |
| <input type="checkbox"/> | <input type="checkbox"/> | Is a instruction manual included? |
| <input type="checkbox"/> | <input type="checkbox"/> | Is there a type label and can it be read easily? |
| <input type="checkbox"/> | <input type="checkbox"/> | Are there all security-related labels and can they be read easily? |

Results of the check of single safety-related components (product-related annex)

- | | | |
|--------------------------|--------------------------|--|
| <input type="checkbox"/> | <input type="checkbox"/> | Has the product related annex been used? |
| <input type="checkbox"/> | <input type="checkbox"/> | Requirements fulfilled? |

Remedy of defects

- | Ja | Nein | |
|--------------------------|--------------------------|--|
| <input type="checkbox"/> | <input type="checkbox"/> | Wurden alle festgestellten Mängel behoben? |
| <input type="checkbox"/> | <input type="checkbox"/> | Have all defects been repaired? |
| <input type="checkbox"/> | <input type="checkbox"/> | Have all defective components been changed? |
| <input type="checkbox"/> | <input type="checkbox"/> | Is the protection equipment working properly? |
| <input type="checkbox"/> | <input type="checkbox"/> | Is there the original type label and can it be read easily? |
| <input type="checkbox"/> | <input type="checkbox"/> | Are there all security-related labels and can they be read easily? |

Changed components

Finish of the test

According to MDD, the safety-related and the hygienic condition of a medical device must be proved or be reconstituted respectively, before it can be used again. Having purified a medical device, which is used conventional, it has to be cleaned, disinfected and if necessary sterilized. Therefore all related working steps (see also RKI-instructions) and tests of proof the technical and functional security have to be performed. Moreover, attention has to be paid to an accurate packaging.

- | Ja | Nein | |
|--------------------------|--------------------------|--|
| <input type="checkbox"/> | <input type="checkbox"/> | Has the device been cleaned? |
| <input type="checkbox"/> | <input type="checkbox"/> | Has the device been disinfected? |
| <input type="checkbox"/> | <input type="checkbox"/> | Has the technical and functional security been reconstituted and proven? |
| <input type="checkbox"/> | <input type="checkbox"/> | Has the device been packed properly? |

Responsible: _____

Date: _____

Signature: _____

Stand: 06/2013