



# VacuAide 7325

## EASY. RELIABLE. EFFICIENT.

The VacuAide 7325 suction device is economical, maintenance free and powerful at the same time. The unit has been designed for all types of respiratory suction therapy for use in the home, nursing homes as well as special clinical applications.

The product range includes two versions: a mains-powered model 7325D and the model 7325P, which is equipped with a long-lasting rechargeable battery.

### Powerful technology in a compact design

By using state-of-the-art components, our engineers have succeeded in developing a lightweight, small footprint, space-saving device with a suction range of 50 to 550 mmHg at a flow rate of 27 l/min. This compact VacuAide 7325P is equipped with the latest generation lithium-ion battery. The battery duration lasts up to 45 minutes depending on therapy usage and only takes 3 hours to fully charge.



### New design and product features

The modern housing design ensures that the collection canister is easy to install and to remove. For convenience and hygiene an integrated tube holder allows the user to place the suction tube after each suction procedure. The unique adjustable vacuum gauge is surrounded in the housing to protect against falling objects and impact. LED control lights are displayed on the top of the device for quickly checking the units charge status.



### TECHNICAL DATA

Dimensions (H x W x D)		17.5 x 17.8 x 18 cm
7325P*	with rechargeable battery	
7325D	mains only operated	
Weight		
Model 7325P	1.72 kg	
Model 7325D	1.54 kg	
Vacuum Range		
50 to 550 mmHg		
Air Flow		
27 l/min		
Battery Running Time (for 7325P)		
up to 45 minutes		
Charging Time (for 7325P)		
3 hours		
Operating Temperature Range		
0° C - 40° C		
Electrical Requirements		
100 - 240 V AC, 50 / 60 Hz; 12 V DC		
Battery Warranty		
1 Year		
Device Warranty		
2 Years		

\*The 7325P-version is supplied with personal carry bag and 12V power adapter

### Features of the new 7325-Series

- Long-life lithium-ion battery with 1 year warranty
- Up to 45 minutes running time
- Fast charging time - only 3 hours!
- Space-saving design, lightweight and powerful
- Tube holder integrated
- Stepless vacuum regulation from 50 to 550 mmHg
- High suction power – up to 27l/min (free flow)
- Functional personal carry bag included (P-version)

### Hygienic suction - protection for the user

A sophisticated safety system prevents liquids from entering the device. The double protection system essentially consists of two components: a filter that closes automatically when it comes into contact with liquid, and a mechanical overflow protection that is located in the lid of the 800ml standard canister.

### Unlimited mobility

In addition to a functional personal carry bag, the standard scope of delivery of the VacuAide 7325P also includes a 12-Volt vehicle auxiliary cable.



### Collection canister – individually Adaptable

The standard contents of delivery with the VacuAide 7325 includes an 800 ml collection canister is intended for same patient use. A reusable autoclavable 1200ml collection canister is available as an option. In a clinical environment, the 7325 unit can accommodate Flovac collection canisters and disposable liners.

### Common VacuAide accessories for efficient Inventory

For greater economic efficiency and simple material planning, all associated collection canisters are compatible with all Drive DeVilbiss VacuAide ranges.



#### Drive DeVilbiss International

DeVilbiss Healthcare GmbH  
Kamenzer Straße 3  
68309 Mannheim, Germany

+44 (0) 121 / 521 3140  
+49 (0) 621 / 1 78 98 0  
kontakt@devilbiss-hc.com

CE 0044  
www.drivedevilbiss-int.com

**drive**  
DeVilbiss  
HEALTHCARE